P1752.2 **Cardiorespiratory Subgroup Teleconference**

March 23, 2023 **UTC 4:00 PM**

IEEE STANDARDS ASSOCIATION

Working Group Sponsored by IEEE Engineering in **Medicine & Biology Standards Committee**



P1752.2 **Cardiorespiratory Subgroup**

Agenda March 23, 2023

Attendance & Introductions Items:

Other business

 \checkmark Preparation (session March 7, 2023): Review JSON schema coding Informed by IEEE 1751.1 schemas Esp familiarize with "Metadata" and "Utility" schemas Consider: *Samples -> Schemas* Assessing dependencies: Organization of CR schema

✓ Interests for schema groupings ✓ Schema timeline / initial draft target ✓ Next meeting : April 27, 2023



Morbidity & mortality ... WHO category of highest prevalence

Preventative medicine ... overall intervention of highest impact

Health economics ... as yet unmet need : semantic interoperability contextuality

Impactful Healthcare Relevance

Cardiovascular disease **#1**

heart failure

cerebrovascular

dysrhythmias

Physical activity & mvt


CV fitness *

Apps

Wearables & ext. detectors

Internal/implant sensors*

PRSteiner - Dartmouth





Morbidity & mortality ... WHO category of highest consequence

Preventative medicine ... overall intervention of highest impact

Health economics ... as yet unmet need : -> semantic interoperability \rightarrow contextuality

Impactful Healthcare Relevance

Cardiovascular disease **#1**

cerebrovascular

heart failure

dysrhythmias

Physical activity & mvt

CV fitness *

Apps

Wearables & ext. detectors


Internal/implant sensors*




IEEE 1752 Standard for Mobile Health Data












Expansion of Mobile Health Data

Standard Patch Rhythm Monitor



- 2 electrodes
- Continuous data up to 2 weeks
- Atrial and ventricular signals
- Morphologic analysis (limited)
- Time of day, clustering
- Dysrhythmia auto-triggering
- Symptom log
- Option for real-time alerts
- Encrypted communications





© 2013 Rhythm Technologies, inc.

bpm, the longest lasting 35.3 secs with an avg rate of 133 bpm. Atrial Fibrillation

triggered/diary event. Isolated SVEs were frequent (5.4%, 76752), SVE Couplets were occasional (3.7%, 26323), and SVE Triplets were occasional (1.7%, 7781). Isolated VEs were rare (0 to <1.0%, 5154), VE Couplets were rare (0 to <1.0%, 19),

Pauses, AV Block, and Atrial Fibrillation were detected within ±45 seconds of patient

and VE Triplets were rare (0 to <1.0%, 1). Ventricular Bigeminy and Trigeminy were

occurred (12% burden), ranging from 45-212 bpm (avg of 84 bpm). 3 Pause(s)

occurred, the longest lasting 5.4 secs (11 bpm). Suprventricular Tachycardia,

present, MD notification criteria for Rapid Atrial Fibrillation and AV Block met -

notified RN on 03/15/2013 at 12:00 pm CST.

S/N: HTEST12345

up to 16 beats rate range of 135-182/min. 8- Atrial fibrillation

to correlate with SVT, sinus pauses, AV Block and atrial


fibrillation.

with a burden of 12% with rates of 45-212/ min, 9- Symptoms appear

Signed by Dr. Example Physician on 03/17/13 at 03:41 PM (CT)

Mobile Health Data: Example Wearable "Long-term" Monitoring

Over 50% of arrhythmias are missed by Holter monitors.²




Over


Detected by Zio XT in 14 days²


7%


Detected by Holter in 2 days²



Mobile Health Data: Wearable "Long-term" Monitoring










IEEE 1752 Standard for Mobile Health Data
















Expansion of Mobile Health Data overlapping into the Digital Biomarker Space...





IEEE 1752 Standard for Mobile Health Data





Expansion of Mobile Health Data overlapping into the Digital Biomarker Space...





Morbidity & mortality ... WHO category of highest consequence

Preventative medicine ... overall intervention of highest impact

Health economics ... as yet unmet need : -> semantic interoperability \rightarrow contextuality

Impactful Healthcare Relevance

Cardiovascular disease **#1**

cerebrovascular

heart failure

dysrhythmias

Physical activity & mvt


CV fitness *

Apps

Wearables & ext. detectors

Internal/implant sensors*





\rightarrow	
	No. N. S. Marine and

Parameter	Value
Volume	®1.4 cm ³
Mass	∲3.4 g
Dimensions H x W x D	∲45.1 mm x ∲8.0 mm x ∲4.2 mm
Surface area of device electrode	
Distance between the electrodes, centroid-to- centroid	∲40 mm

